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
**Caution:** Clean the rubber inside the oximeter and the finger to be inserted with surgical spirit before using the device. The rubber does not contain toxins and is well tolerated by human skin.

- Put the finger between the rubber-coated halves of the device and close the clip carefully.
- Press the power button on the front of the device once.
- Do not move the finger during the measurement and avoid any other movements.
- Wait until the device has captured the signal, then read the measured values from the display.

### Display

Switch on the oximeter. After capturing the measured values, the pulse rate, SpO<sub>2</sub> saturation and a pulse graph will appear on the LED display.

### Product features

- If no signal is measured for more than 5 seconds, the device will switch off automatically.

### Maintenance and storage

- Clean the surface of the device before each examination.
- Remove the batteries if the device is taken out of service for a longer time.
- Adhere to the national and local legislation when handling and/or disposing of batteries.

### Accessories

- Carrying strap
- Instructions for use

### Troubleshooting

#### Problem

Oxyhaemoglobin and pulse are not displayed.

#### Causes

1. The finger was not inserted correctly.
2. The oxygen saturation of the patient is below the measuring range.
3. Nail varnish may have been applied.

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Lagerungstemperatur: -10-40°C

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### Solutions

1. Insert finger again.
2. Try again several times, otherwise go to hospital for an accurate measurement.
3. Remove nail varnish.

#### Problem

Oxyhaemoglobin and pulse are displayed irregularly.

#### Causes

1. The finger was not inserted deeply enough.
2. The finger was trembling, or the patient moved.

#### Solutions

1. Insert finger again.
2. Avoid movements.

#### Problem

Oximeter does not switch on.

#### Causes

1. The batteries are discharged.
2. The batteries were not inserted correctly.
3. The oximeter might be damaged.

#### Solutions

1. Replace batteries.
2. Take batteries out and insert them again.
3. Contact your dealer.

#### Problem

The display switches off suddenly.

#### Causes

1. The device switches off automatically after 5 seconds if no signal is received.
2. The battery voltage is no longer sufficient.

#### Solutions

1. This is normal.
2. Replace batteries.

### Technical data

Display: LED

SPO<sub>2</sub>: Measuring range: 70-99%,

Deviation: ±2%

Pulse: Measuring range: 30-235 BPM,

deviation: ±2 BPM or ±2%

Pulse intensity: Bar graph

Power display: 2 AAA Alkaline,

Battery life: up to 30 hours

Ambient conditions: Temperature: 5-40°C,

Storage temperature: -10-40°C

Relative humidity during operation: 15%-80%

Relative humidity during storage: 10%-80%

### Precautions:

- Do not use the oximeter near CT or MRI devices.
- Danger of explosion! Do not use the oximeter in an environment where there is a risk of explosion.
- The oximeter is designed as a utility to complement the assessment of a patient. It has to be used in conjunction

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with other methods of assessment of clinical parameters and symptoms.

- Regularly check whether the oximeter fits correctly.
- Read these instructions carefully before using the device.
- The oximeter is not provided with an alarm function. It is therefore not suitable for the continuous monitoring of the oxygen saturation of blood.
- The function of the oximeter and the correct fit of the sensor should be checked at least every 4 hours.
- Autoclaving (i. e. steam pressure sterilization), ethylene oxide sterilization, or immersing the device in liquid can lead to wrong measuring values.
- Certain functional failures of haemoglobin may affect the values.
- Certain dyes used for medical diagnosis may affect the measured values.
- Bright light may also affect the values. If necessary, shield the device with a cloth or a similar object.
- Intense movements of the patient may distort the measuring results.
- High frequencies may distort the measuring results.
- Venous pulsation may also distort the measuring results.
- Placing the sensor near a pressure bandage, an arterial or venous access may affect the values.
- Hypotension, vasoconstriction, anaemia, hypothermia, cardiac arrest, and shock may affect the values.
- Nail varnish and false nails can also distort the values.

### 2 year's warranty

Batteries, rechargeable batteries and electrical devices must not be disposed of with the normal household waste! All consumers are legally obliged to submit all batteries, rechargeable batteries and electrical devices, whether they contain pollutants or not, to a collection point in their community or district or to return them to a dealer, so that they may be disposed of an environmentally compatible way. Only submit batteries and rechargeable batteries in a discharged state!

Beijing Choice Electronic Technology Co., Ltd.  
Room 4104, No. A12 Yuqun Road, Haidian District  
100143 Beijing  
PEOPLE'S REPUBLIC OF CHINA

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## GB Instructions for use

### General description

The Fingertip Pulse Oximeter is designed to measure the oxygen saturation of arterial haemoglobin (SpO<sub>2</sub>%) and the pulse rate.


**Caution:** The Fingertip Pulse Oximeter is intended as an additional utility for assessing patients. The displayed values have to be considered in conjunction with other clinical parameters and symptoms.

### Fields of application

The Fingertip Oximeter can be used to measure the haemoglobin saturation and heart rate in humans by using the finger. The product is suitable for use in hospitals (including clinical use in the fields of internal medicine, surgery, anaesthesia, paediatric care, intensive care, etc.), oxygen bars, and sociomedical organizations as well as for medical care in sports.

**Caution:** The device can be used before or after training. A use during training is not recommended.

### Operational diagram

1. Light and infrared transmitter
2. Light and infrared receiver

### Operation

#### Insert the Batteries

- Insert the two batteries with correct polarity.
- Close the lid of the battery compartment as shown below.

**Caution:** Pay attention to the correct polarity of the batteries – otherwise the device may be damaged. Never insert the batteries forcefully. Remove the batteries from the device if it is taken out of service for a longer time.  
- Open the clip by pressing it as shown in the figure.

Insert a finger with the fingernail showing upwards.

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